

## RED HM TAQ DNA POLYMERASE

**Source:** An *E. Coli* strain that carries HM Taq DNA polymerase gene.

**Application:** Routine PCR amplification of DNA fragments upto 6 kb from genomic DNA.

**Storage Buffer:** 12.5 mM Tris-HCl (pH 8.5) with optimized concentration of DTT, EDTA and 50%(V/V) Glycerol.

**Shipping Condition:** Shipped in ice packs.


**Storage conditions:** Storage at -20°C is recommended.

### Package Contents:

- ✓ HM Red Taq DNA Polymerase
- ✓ 10×R-Taq Buffer
- ✓ 20mM MgCl<sub>2</sub>

**10×R-Taq buffer composition:** 200 mM Tris-HCl (pH 9.2) with optimized concentration of KCl & MgCl<sub>2</sub> along with NP-40.

- **Product Code:** RT002
- **Quantity;** 100 unit
- **Lot No.:**
- **Expiry:**




### PCR Reaction

Components	Volume(for 20µl reaction)
10 × R –Taq Buffer	2 µl
10 mM dNTPs	0.6 µl
10 µM Forward Primer	0.4 µl
10 µM Reverse Primer	0.4 µl
HM Red Taq DNA Polymerase	1 µl
DNA Template	Variable
Nuclease Free Water	14.6 µl

### PCR Program(General)

Step	Temperature	Time
Initial denaturation	95°C	1-5 min
25-35 cycles (Annealing)	95°C	30 sec
	45-68°C	1-2 min
	72°C	30 sec- 2 min
Final Extension	72°C	1-5 min
Hold	4°C	∞



### Performance Test :

Quality control analysis of Red Taq DNA polymerase (Lot No.: PT13/243): Red Taq used for PCR of plant genomic & bacterial plasmid DNA template.

### Ordering Information:

GeNext Genomics Pvt.Ltd.  
103, Abhyankar Nagar,  
Nagpur  
Ph.No.08888803973

**Note:** No loading buffers or tracking dyes required. Samples may be added directly to an agarose gel after PCR and visualized.



### Features:

#### High Purity & specificity

Direct loading of PCR products without the need of loading dye for visualization